# Master all your challenges

*initial*<sup>TM</sup>

Press-over-Metal  
Press-over-Zircon

- EN Technical Manual
- DE Technische Arbeitsanleitung
- FR Manuel Technique
- IT Manuale Tecnico
- ES Manual Técnico

JULY 2021

**GC**

EN

DE

## Table of Contents

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Page

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EN GC Initial IQ, One Body, Press-over-Metal System Colour Chart

DE GC Initial IQ, One Body, Press-over-Metal System Farbkarte

Description	No. of Bottles	Colour Code
Effect Dentin	7	DA-1
		DA-2
		DB-1
		DB-2
		DC-1
		DC-2
		DD-2
Effect Bleach Dentin	3	A0
		B0
		B00

EN GC Initial IQ, One Body, Press-over-Zircon System Colour Chart

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EN GC Initial IQ • One-Body,  
Press-over-Metal Opaque Layering

DE GC Initial IQ • One-Body,  
Press-over-Metal Opaquerschichtung


EN Cast or CAD/CAM made metal frameworks are finished using tungsten carbide metal burs or ceramic-bonded grinding instruments and oxidized according to alloy manufacturer's instructions.

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EN Apply the first Opaque layer (wash) with the normal GC Initial MC Paste Opaques in a thin layer using a flat brush. Before firing fluorescent crystals are applied on the Opaque surface.

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EN After firing these Fluo Crystals are the base for a very important equal layer of the 2nd Opaque firing.

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EN \*Advised temperature when firing on NPA (Non Precious Alloys)

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EN GC Initial IQ • One-Body, Press-over-Zircon  
Power Frame Modifiers

DE GC Initial IQ • One-Body, Press-over-Zircon  
Power Frame Modifiers

EN For optimized shade reproduction, it is recommendable to use a pre-coloured Zirconium dioxide framework.

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EN The pre-coloured Zirconium frame structures can be changed in colour and fluorescence by applying a thin wash layer of Power Frame Modifiers.  
Remark: Apply the Power Frame Modifier in a thin layer.

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EN It is unfortunately impossible to offer a general colour combination chart for Power Frame Modifiers with the 16 Vitapan Classic® shades as the existing Zr frame colouring liquids do vary too much from manufacturer to manufacturer.

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EN GC Initial IQ • One-Body, Press-over-Zircon  
Power Frame Modifiers

DE GC Initial IQ • One-Body, Press-over-Zircon  
Power Frame Modifiers

EN GC Initial IQ Power Frame Modifiers in the posterior regions. Aesthetical critical areas should be supported by extra colour from the Power Frame Modifiers.

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


- EN Sprueing of the modellation:  
The wax sprues (3.0-3.5mmØ) should be app. 5-6 mm in length and should not taper towards the pattern (Rounded out – no sharp edges, consider the direction of flow of the ceramic material: same flow direction).  
The patterns are waxed carefully onto the ring base former in an angle of 45° (Distance between the individual objects: min. 3mm).  
Remark: Multi-unit bridge constructions should be made with a crossbar instead of direct sprueing.

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EN Wax weight: Please weigh the objects before investing. To get the right wax-weight, the frame structure has to be weighed before and after the anatomical wax modellation.

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
EN The investing is carried out with GC MultiPressVest®, specially developed, fine grained phosphate bonded investment (For details please refer to the special instruction for use).

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- EN Mix GC MultiPressVest® according the instructions for use:  
premix powder/liquid by hand with a spatula until a uniform mixture.

Mixing ratio:

Small ring size: 100g powder – 22ml liquid

Medium ring size: 200g powder – 44ml liquid

Large ring size: 300g powder – 66ml liquid

Mix for 60 sec. under vacuum (320 – 420 rpm)

After filling the investment ring carefully, put on the mould base former turning it gently (Note: the investment needs to extrude easily through the opening). After setting remove the mould base by turning carefully. Push the investment cylinder out of the silicone ring. Smooth the underside by using a plaster knife (90° angle, stable position).

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# EN Investing

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## EN Notes:


- Check the correct fit of the silicone cylinder before investing.
- The base of the cylinder must be absolutely flat (90° angle), thus the mould stands perfectly in an upright, stable position in the injection furnace.
- Investment residues may not enter the sprues. Clean them carefully before preheating.
- Investment materials contain quartz powder. Avoid the inhalation of dust.
- Please observe the corresponding Instruction for Use for GC MultiPressVest®.

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### EN Preheating:

#### A.: Quick heating method:

After 20 min. put the investment (the funnel pointing downwards) centrally in the furnace, pre-heated to 850°C (Pre-heating: 60min).

If there is more than 1 ring in the furnace, each holding time should be extended by 10 min per add. cylinder.

- GC Initial IQ Press Pellets should not be preheated.
- Alox-press plungers can be put straight at 850°C in the oven.
- One way press plungers do not need to be pre-heated.
- All Alox-press plungers need to be clean before pressing (Use glassbeads for sandblasting).
- Attention: short timing and distance between preheating furnace and pressing is needed!

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EN All data are general references of in house tests of calibrated press furnaces. Please calibrate your press furnaces regularly and test the final temperature. In case of too low pressing temperatures, the CTE-value of the press ceramic may change, resulting in cracking and chipping.

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EN After cooling, please mark the length of the plunger on the investment ring and cut separate along the length of the mark using a suitable cut-off disc (sintered diamond disc).

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EN Afterwards carefully break apart the 2 sections.

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EN Roughly sandblast the injected patterns with glass beads (4 bar, 50µm), after this clean the objects carefully and completely (2 bar, 50µm).

Attention: Do not use Aluminium oxide!

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## EN Devesting and Finishing

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EN Cut of the sprues carefully using a suitable cut-off disc (fine diamond disc, no pressure and low rotation). The ceramic objects are prepared according to the technique used (fine diamonds, low RPM, avoid heating up to prevent cracking). The contact areas metal frame / ceramic (e.g. ceramic shoulder) and the Zircon frame structure / ceramic have to be treated with caution. In case of small failures, a correction powder of each ingot colour is available in the sets.

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
## EN Lustre Pastes NF Application - Basic Step by Step Instruction

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**EN** Important Information: Please always stir up material well in the jars before use. Before applying the Lustre Coating, the ceramic restoration is shaped and contoured using adjusted diamond burs. After prepolishing with a ceramic polisher the finished crowns & bridges are sand-blasted with 50 microns Alu-Oxid at 1,5 bar pressure.

The Lustre Paste Neutral is coated on the entire surface of the sandblasted restoration. This coating is clearly thicker than what we know from a normal glaze firing.

Note: These Lustre Pastes may not come in contact with water.

The requested V-Classic Shade is created by a more or less intensive application of the Lustre Pastes "Body A - D". This process, which only takes 2 minutes, is everything which is necessary to create V-Classic shades.

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



EN To build up a high aesthetic solution for the patient use our advanced step by step instruction:

In order to achieve the higher grey value of the incisal area with darker teeth (A3.5; A4; C3; C4), 2 grey Effect Shades can be used: L-3 (Light Grey) and L-4 (Dark Grey). For eventual further individualisation, you can use the whitish Lustre Pastes NF: L-1 (Vanilla) and L-2 (White) and L-V (Value) for brightening the restorations. 2 bluish Effect Lustre Pastes NF: L-5 (Light Blue) and L-6 (Dark Blue) for a bluish opalescent appearance in the incisal area.

Greenish Effect Lustre Paste NF: L-8 (Olive) for a greenish effect. Can be mixed with all Lustre Pastes NF.

In order to achieve an even higher „Transparency“ the 2 greyish Effect Lustre Pastes NF L-3/L-4 can be used as well. All GC Initial Spectrum Stains / GC Initial MC INvivo Stains can be used to achieve even more possibilities of individualization.

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EN The requested surface texture / smoothness of the Lustre Paste layer can be obtained by soft vibrating or condensing the restoration.

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EN No matter if you choose the basic way or the advanced way for high aesthetic results, the steps to finalize your work are still the same.

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
EN GC Initial Press-over-Metal is a veneering ceramic for all types of precious and non-precious alloys. Make sure that the CTE of the alloy you use for casting or milling the substructure stays strictly within the range of 13,8 - 14,9 · 10<sup>-6</sup> K<sup>-1</sup> (25-500 °C). Calibrate your furnace according to the manufacturer's instruction in order to obtain the best results and to be able to make optimal use of the properties. As each dental ceramic is sensitive to contamination, take care that you work in a clean environment.

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
FR GC Initial IQ, One Body, Press-over-Metal System  
Tableau des teintes

IT GC Initial IQ, One Body, Press-over-Metal System  
Tabelli Colori


Description	No. of Bottles	Colour Code
Effect Dentin	7	DA-1
		DA-2
		DB-1
		DB-2
		DC-1
		DC-2
		DD-2
Effect Bleach Dentin	3	A0
		B0
		B00

FR GC Initial IQ, One Body, Press-over-Zircon System  
Tableau des teintes

IT GC Initial IQ, One Body, Press-over-Zircon System  
Tabelli Colori

Description	No. of Bottles	Colour Code
Effect Dentin	7	DA-1
		DA-2
		DB-1
		DB-2
		DC-1
		DC-2
		DD-2
Effect Bleach Dentin	3	A0
		B0
		B00

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
FR GC Initial IQ • One-Body,  
Press-over-Metal - couche opaque

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FR GC Initial IQ • One-Body, Press-over-Zircon, Power Frame Modifiers

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**FR** GC Initial IQ • One-Body, Press-over-Zircon, Power Frame Modifiers

**IT** GC Initial IQ • Una sola massa, Press-over-Zircon, i Power Frame Modifier

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**EN Medical Device Regulation Requirement**

Undesired effects - Reporting: If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: [https://ec.europa.eu/growth/sectors/medical-devices/contacts\\_en](https://ec.europa.eu/growth/sectors/medical-devices/contacts_en) as well as to our internal vigilance system: [vigilance@gc.dental](mailto:vigilance@gc.dental)  
In this way you will contribute to improve the safety of this product.

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RESPONSIBLE MANUFACTURER IN CANADA  
GC AMERICA INC.  
3737 W. 127th Street, Alsip, IL 60803 U.S.A.

GC AMERICA INC.  
3737 West 127th Street, Alsip, IL 60803 U.S.A.  
Tel: +1-708-597-0900  
[www.gcamerica.com](http://www.gcamerica.com)

GC Germany GmbH  
Seifgrundstraße 2, D-61348 Bad Homburg  
Tel. +49.61.72.99.59.60, Fax. +49.61.72.99.59.66.6  
info.germany@gc.dental, <https://europe.gc.dental/de-DE>

GC ITALIA S.r.l.  
Via Calabria 1, I-20098 San Giuliano Milanese  
Tel. +39.02.98.28.20.68, Fax. +39.02.98.28.21.00  
info.italy@gc.dental, <https://europe.gc.dental/it-IT>

GC UNITED KINGDOM Ltd.  
Coopers Court, Coopers Court, Newport Pagnell, UK-Bucks. MK16 8JS  
Tel. +44.1908.218.999, Fax. +44.1908.218.900  
info.uk@gc.dental, <https://europe.gc.dental/en-GB>

GC FRANCE s.a.s.  
8 rue Benjamin Franklin, F-94370 Sucy en Brie Cedex  
Tel. +33.1.49.80.37.91, Fax. +33.1.45.76.32.68  
info.france@gc.dental, <https://europe.gc.dental/fr-FR>

GC IBÉRICA - Dental Products, S.L.  
Edificio Codesa 2  
Playa de las Américas, 2, 1º, Of. 4, ES-28290 Las Rozas, Madrid  
Tel. +34.916.364.340, Fax. +34.916.364.341  
comercial.spain@gc.dental, <https://europe.gc.dental/es-ES>

GC AUSTRIA GmbH  
Tallak 124, A-8103 Gratwein-Strassengel  
Tel. +43.3124.54020, Fax. +43.3124.54020.40  
info.austria@gc.dental, <https://europe.gc.dental/de-AT>

GC Europe NV - Benelux Sales Department  
Researchpark Haasrode-Leuven 1240  
Interleuvenlaan 33, B-3001 Leuven  
Tel. +32.16 74.18.60  
info.benelux@gc.dental, <https://europe.gc.dental/fr-BE>

GC EUROPE N.V. - East European Office  
Siget 19b, HR-10020 Zagreb  
Tel. +385.1.46.78.474, Fax. +385.1.46.78.473  
info.eeo@gc.dental, [www.eeo.gceurope.com](http://www.eeo.gceurope.com)

GC NORDIC AB - Finnish Branch  
Lemminkäisenkatu 46, FIN-20520 Turku  
Tel. +358.40.900.07.57  
info.finland@gc.dental, <https://europe.gc.dental/fi-FI>

GC NORDIC AB - GC Nordic Danish Branch  
Scandinavian Trade Building  
Gydevang 34-41, DK-3450 Allerød  
Tel. +45 51 15 03 82  
info.denmark@gc.dental, <https://europe.gc.dental/da-DK>

GC NORDIC AB  
c/o Lundin Revisionbyrå  
Erik Dahlbergsgatan 11B, SE-411 26 Göteborg  
Tel. +46.768.54.43.50  
info.nordic@gc.dental, <https://europe.gc.dental/sv-SE>

GC AUSTRIA GmbH - Swiss Office  
Zürichstrasse 31, CH-6004 Luzern  
Tel. +41.41.520.01.78, Fax. +41.41.520.01.77  
info.switzerland@gc.dental, <https://europe.gc.dental/de-CH>

GC AUSTRALASIA DENTAL PTY LTD  
1753 Botany Rd, Banksmeadow NSW 2019, Australia  
Tel: +61-2-9301 8200, Fax: +61-2-9316 4196

GC SOUTH AMERICA  
Rua Heliadora, 399, Santana - São Paulo, SP, BRASIL  
CEP: 02022-051 - TEL: +55-11-2925-0965 - CNPJ: 08.279.999/0001-61  
RESP. TÊC: Erick de Lima - CRO/SP 100.866

GC ASIA DENTAL PTE. LTD.  
5 Tampines Central 1, #06-01 Tampines Plaza 2, Singapore 529541  
Tel: +65 6546 7588

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EN V-Shade is referring to the Vita® Shade Guide. Vita® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.