everStick®POST

Individually formable glass fibre root canal posts

DISTRIBUTED BY

GC CORPORATION

76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

GC AMERICA INC.

3737 West 127th Street, Alsip, IL 60803 U.S.A.

TEL: +1-708-597-0900

GC ASIA DENTAL PTE. LTD.

11 Tampines Concourse, #03-05, Singapore 528729

TEL: +65 6546 7588

GC AUSTRALASIA DENTAL PTY. LTD

1753 Botany Rd, Banksmeadow, NSW 2019, Australia

TEL: +61 2 9316 4499

MADE IN FINLAND

Fibre type: Silanated E-glass fibre impregnated with bis-GMA and PMMA

Form: Unidirectional fibre bundle

Diameter: ~0.9 mm; ~1.2 mm; ~1.5 mm

Caution: US Federal law restricts this device to sale by or on the order of a dentist.
U7204 - 951332

everStick®POST

- EN INDIVIDUALLY FORMABLE GLASS FIBRE ROOT CANAL POSTS
- DE INDIVIDUELL FORMBARE GLASFASERN FÜR WURZELKANALSTIFTE
- FR TENON RADICULAIRE ANATOMIQUE EN FIBRE DE VERRE
- IT PERNI ENDOCANALARI IN FIBRA MODELLABILI INDIVIDUALMENTE
- ES POSTES DE FIBRA DE CANAL RADICULAR MODELADOS INDIVIDUALMENTE
- NL INDIVIDUEEL TE VORMEN GLASVEZELWORTELKANAALSTIFTEN
- DA INDIVIDUELLE OG AVANCEREDE RODSTIFTER
- SV INDIVIDUELLT ANPASSADE ROTKANALSTIFT I GLASFIBER
- PT PINOS DE FIBRA DE VIDRO MALLEÁVEIS PARA CONDUTOS RADICULARES
- EL ΕΝΔΟΡΡΙΖΙΚΟΙ ΑΞΟΝΕΣ ΑΠΟ ΙΝΕΣ ΥΑΛΟΝΗΜΑΤΩΝ ΜΕ ΕΞΑΤΟΜΙΚΕΥΜΕΝΗ ΔΙΑΜΟΡΦΩΣΗ
- FI JUURIKANAVANASTA

GC

EN	Individually formable glass fibre root canal posts	5	DA	Individuelle og avancerede rodstifter	41
DE	Individuell formbare Glasfasern für Wurzelkanalstifte	11	SV	Individuellt anpassade rotkanalstift i glasfiber	47
FR	Tenon radiculaire anatomique en fibre de verre	17	PT	Pinos de fibra de vidro maleáveis para condutos radiculares	53
IT	Perni endocanalari in fibra modellabili individualmente	23	EL	ΕΝΔΟΡΡΙΖΙΚΟΙ ΑΞΟΝΕΣ ΑΠΟ ΙΝΕΣ ΥΑΛΟΝΗΜΑΤΩΝ ΜΕ ΕΞΑΤΟΜΙΚΕΥΜΕΝΗ ΔΙΑΜΟΡΦΩΣΗ	59
ES	Postes de fibra de canal radicular modelados individualmente	29	FI	Juurikanavanasta	65
NL	Individueel te vormen glasvezelwortelkanaalstiften	35			

Prior to use, carefully read the instructions for use.

INSTRUCTIONS FOR USE:

The everStickPOST is an adaptable, polymer (PMMA) and resin-impregnated (bis-GMA) unpolymerised glass fibre post. Polymerising this material produces a post with high flexural strength and elasticity very similar to the natural elasticity of dentine. Consequently, the stress of occlusion will be evenly distributed on the root structure. Adhesive and micromechanical bonding to both resin cement and composite ensures a strong bond to the root canal and the composite core.

When everStickPOST's are used, preparation of the root canals need not be as extensive as with traditional posts. Thus, the dentine can be saved and the risk of perforation is reduced because the canal preparation is

minimised. The pulp chamber of the root canal can be completely filled with fibres instead of cement. When the post is adapted to the morphology of the canal and the root canal is filled with fibres, the adhesive surface and the strength of the most critical part of the tooth are maximised.

The unique properties of everStickPOST glass fibre posts also make it possible to use the posts in curved and oval root canals as well as in very large canals, where several posts of different lengths and diameters can be placed in the same canal. Also, everStickPOST's can be used in traditionally prepared and enlarged root canals.

CONTRA-INDICATION

In rare cases the product may cause sensitivity in some people. If any such

reactions are experienced, discontinue the use of the product and refer to a physician.

CLINICAL PROCEDURES

Prior to the use of everStickPOST, the root canal(s) of the tooth should be endodontically treated and filled in accordance with generally approved methods. In order to make the best use of the properties of everStickPOST, it is recommended to apply tissue saving principles during the preparation of the root canal.

SELECTION OF THE everStickPOST SIZE

The everStickPOST fibres are available in three diameters: 0.9, 1.2, and 1.5 mm. The most suitable size can be selected for canals of many different sizes and shapes. In large canals and root canal openings, use of two or more posts is recommended.

ROOT CANAL PREPARATION

1. Remove 2/3 of the length of the root canal filling material or at least the height of the clinical crown from the canal. For example, gutta-percha can be removed with a suitable size Gates Glidden bur without enlargement of the canal. All of the root canal filling material should be removed from the length of the preparation. Leave at least 3–5 mm of gutta-percha at the apex of the root. Rinse the canal with water, and dry carefully with paper points. The working area should be isolated from moisture as well as possible. Use of rubber dam isolation is highly recommended.
2. Measure the depth of the prepared canal, using, for example, an endodontic instrument or a periodontal probe. Also estimate the height of the coronal structure required.

3. Open the everStickPOST foil bag. Cut the required number of posts from the silicone strip, using scissors. Close the foil bag with its sticker, and place the closed bag in the refrigerator.

PLACEMENT PROCEDURE of everStick-POST

4. Mark the measured length of the post on the protective paper. Pre-cut the post together with the silicone to a suitable length. Use sharp scissors.
5. Use tweezers to remove the post from the silicone. Check the length and suitability of the everStickPOST by inserting it into the root canal. Always use tweezers to handle the post.

CLINICAL TIP: Dipping the tweezers in a

drop of light curing resin (e.g., StickRESIN) prevents them from sticking to the everStick-POST. The resin will also enhance the bonding between the main post and the additional post(s) attached to it. Avoid letting resin flow into the root canal.

6. If the post does not reach the necessary depth, taper the end of the post with sharp scissors.
7. Fit the post inside the root canal again. At this point, if necessary, you can shorten the coronal section of the post to an appropriate length with sharp scissors.
8. In the upper portion of an oval or very large canal, it is recommended to use more than one post in order to strengthen the post in areas of greater load. All

additional posts are shaped and attached tightly to the main post both coronally and inside the root canal by means of lateral condensation.

Important: Remove the post from the canal and protect it from the light before cementing.

CEMENTING everStickPOST

To cement the everStickPOST in place, use low-viscosity dual-curing cement. Pay careful attention to the manufacturer's instructions.

CLINICAL TIP: It is important to select dual-curing composite resin cement with low viscosity. Use of highly viscous cement may prevent the unpolymerised post from reaching full depth inside the canal.

9. Follow the instructions of the cement manufacturer to prepare the root canal prior to cementing. Fill the canal with cement, using an intraoral tip. It is important to start filling the canal from the apical region and proceed slowly by moving the syringe steadily upwards until the canal is filled.

NOTE: Do not use a Lentulo spiral to apply the cement – it accelerates the polymerisation process of the composite cements.

NOTE: Covering the post with cement instead of filling the canal prior to inserting the post into the canal may cause air voids and shredding of the individual posts from the post bundle.

10. Slowly insert the post into the canal. You can shape and bend the coronal part of

the post while it is still soft. Be careful not to lift the post at this point. You can remove any excess cement now.

11. Light-cure the post and the cement from above the post, perpendicular to the fibres, for at least 40 seconds.
12. When the post and the cement have cured, you can continue building the coronal part of the tooth, using the preferred method and the composite material best suited to the purpose.

ALTERNATIVE INSTRUCTIONS FOR USE

These instructions can be followed if it is expected that for certain reasons (e.g., highly viscous cement or a long and narrow root canal) the post may not reach the prepared depth during cementing.

First follow steps 1–7 of the above instructions for use; then apply steps A–F below.

- A. In the upper portion of an oval or very large canal, you can use more than one post to strengthen the post in areas of greater load. Any additional posts are shaped and bonded to the main post coronally with a thin layer of light-curing enamel bonding agent* (for example, StickRESIN). Avoid letting enamel resin flow into the root canal. Light-cure it for 20 seconds inside the root canal before removing the post from the canal.
- B. Important: Remove the post from the canal and light-cure it for a total of 40 seconds, making sure that the post is cured on all sides. Fit the post in the canal again, ensuring that it reaches the appropriate depth.

Tip: Locking tweezers will help keep the post in the correct position when it is cemented. If there are other canals in the tooth, prepare posts for them in the same manner.

CEMENTING everStickPOST

- C. Activate the surface of the post(s) thoroughly, using a light-curing enamel bonding agent*, such as StickRESIN. Place the post(s) under a light shield for 3–5 minutes to prevent premature curing. During activation, prepare the canal(s) for cementing as described in the cement manufacturer's instructions. Prior to cementing, thin the resin layer by carefully blowing dry, oil-free air over the surface of the post. Light-cure the post thoroughly for 10 seconds. Use chemically curing or dual-curing cement, paying careful attention to the manufacturer's instructions.
- D. Fill the canal with cement as instructed by the manufacturer. It is important that you start filling the canal from the apical region and proceed slowly by moving the syringe steadily upwards until the canal is filled.
- E. Slowly insert the post into the canal. Add any other necessary posts to the canal(s) in the same manner. Hold the posts in place until the composite cement is sufficiently cured. Light-cure the dual-curing cement according to the composite cement manufacturer's instructions.
- F. After the cement is cured, you can continue building the coronal part of the tooth, using the preferred method and the composite material most suitable for the purpose.

* **NOTE:** The enamel bonding agent used for attaching additional posts and activating the surface of the post must be monomer-based, and it must not contain solvents (acetone, alcohol, water, etc.). The bonding agents in composite cementing kits are not necessarily suitable, because they may contain solvents.

REMOVAL OF AN everStickPOST

The procedures for removal of traditional glass fibre posts can also be applied to the removal of everStickPOST.

STORING: All everStick products should always be stored in a refrigerator (2-8°C; 35-46°F). In addition, you should protect the products from light by keeping them in their sealed foil packages after they have been opened. Higher temperatures and exposure to bright light may shorten the life of everStick

products. Prior to application, the products are taken out of the refrigerator and the foil package opened – but still kept away from bright sunlight or artificial light. When you are cutting the fibre bundle, the rest of the fibre bundle should remain inside the foil package for protection against light. Immediately after a sufficient length is cut for the fibre construction, the foil package is carefully resealed and returned to the refrigerator.

(Shelf life : 2 years from date of manufacture)

PACKAGES

everStick POST INTRO:

5x 2cm post Ø 0.9; 5x 2 cm posts Ø 1.2; 5x 2cm posts Ø 1.5; 5ml StickRESIN bottle

Refills

10 x 2 cm posts; available in size 0.9; 1.2; 1.5

NOTE: These products should be used clinically with care, and the patient should be warned not to abrade the fitting surface, so as to avoid exposing irritation-causing fibres. After the final light-curing of 40 seconds, polymerisation of the fibres will continue for the next 24 hours before they achieve their final strength. Clean the StickSTEPPER and StickCARRIER hand instruments and the StickREFIX D and StickREFIX L silicone instruments by sterilising them before use.

WARNING: Unpolymerised resin can cause skin sensitisation to acrylates in some people. If your skin comes in contact with resin, wash it thoroughly with soap and water. Avoid contact of uncured material with skin, mucous membrane, or eyes. Unpolymerised everStick products may have a slight irritating effect and lead to sensitization to meth-

acrylates in rare cases. The use of powder-free gloves is recommended with everStick products. Polymerize everStick before waste disposal.

US Federal law restricts this device to sale by or on the order of a dentist.

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
Keep away
from sunlight

Temperature
limit







Caution: US Federal law restricts this device to sale by or on the order of a dentist.
U7302 - 951331

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CEMENTING everStickPOST

- C. Activate the surface of the post(s) thoroughly, using a light-curing enamel bonding agent*, such as StickRESIN. Place the post(s) under a light shield for 3–5 minutes to prevent premature curing. During activation, prepare the canal(s) for cementing as described in the cement manufacturer's instructions. Prior to cementing, thin the resin layer by carefully blowing dry, oil-free air over the surface of the post. Light-cure the post thoroughly for 10 seconds. Use chemically curing or dual-curing cement, paying careful attention to the manufacturer's instructions.

- D. Fill the canal with cement as instructed by the manufacturer. It is important that you start filling the canal from the apical region and proceed slowly by moving the syringe steadily upwards until the canal is filled.
- E. Slowly insert the post into the canal. Add any other necessary posts to the canal(s) in the same manner. Hold the posts in place until the composite cement is sufficiently cured. Light-cure the dual-curing cement according to the composite cement manufacturer's instructions.
- F. After the cement is cured, you can continue building the coronal part of the tooth, using the preferred method and the composite material most suitable for the purpose.

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STORING: All everStick products should always be stored in a refrigerator (2-8°C, 35-46°F). In addition, you should protect the products from light by keeping them in their sealed foil packages after they have been opened. Higher temperatures and exposure to bright light may shorten the life of everStick

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