EN

**G-CEM ONE™ Paste Pak**  
SELF-ADHESIVE RESIN CEMENT  
**G-CEM ONE™ ADHESIVE ENHANCING PRIMER**  
TOOTH PRIMER FOR G-CEM ONE

For use only by a dental professional in the recommended indications.

**INDICATIONS FOR USE**  
1. Cementation of all types of all ceramic, resin and metal-based inlays, onlays, crowns and bridges.  
2. Cementation of metal, ceramic, fiber posts, cast post and cores.  
3. Cementation of all ceramic and composite veneers.  
4. Final cementation of crowns and bridges on implant abutments.

**CONTRAINDICATIONS**  
1. Direct pulp capping  
2. Amalgam in contact in patients with known allergies to methacrylate monomer, methacrylate polymer or ethanol.

**IDENTIFICATION OF PARTS (Fig. 1)**  
A. Paste Pak Cartridge  
1. Cover 2. Inner plunger 3. Cartridge tip  
B. Piston release  
4. Piston 5. Adjusting guide 6. Lever 7. Sliding block 8. Piston release

**DIRECTIONS FOR USE**

Mixing Time (sec)	10"
Working Time (23°C / 73°F) (min., sec.) (from start of mixing)	2'30"
Finishing Time (37°C / 99°F, after excess cement removal) (min., sec.)	4'00"

**PASTE PAK CARTRIDGE LOADING**

1. Make sure that the piston is completely retracted into the dispenser (Fig. 2).  
NOTE:  
To retract the piston, push the piston release and at the same time pull the sliding block (on the underside of the dispenser) with fingers.  
2. Load the adhesive resin into the dispenser by moving the cartridge fully forward. Ensure the ▲ mark on the cartridge lines up with the ▲ mark on the Paste Pak Dispenser as illustrated. Then turn the cartridge fully to the right position.  
NOTE:  
To remove the cartridge, retract the piston fully and turn the cartridge back to the marked position.  
3. After cartridge loading, push the piston release forward until it stops (Fig. 4).  
a) Push the piston release forward.  
b) Do not push the upper side.

**A. CEMENTATION TECHNIQUE for all inlays, onlays, crowns, bridges and veneers**

1. TOOTH PREPARATION  
a) Prepare tooth in usual manner.  
b) Clean the prepared tooth with pumice and water. A clean surface is paramount for an optimal adhesion.  
c) Rinse thoroughly with water. Dry thoroughly by gently blowing with oil free air.  
d) Prepared surface should be dry.  
e) In case of non-retentive preparation where extra adhesion is needed, apply G-CEM ONE ADHESIVE ENHANCING PRIMER to the preparation.  
f) In case of veneer or onlay cementation, enamel etching with 35-40% phosphoric acid should be applied for 10-15 seconds and rinsed off, followed by the application of G-CEM ONE ADHESIVE ENHANCING PRIMER to adhesive surface.  
NOTE:  
1. Do not use substances like H<sub>2</sub>O<sub>2</sub>, EDTA and/or disinfectants as these may lower the bond strength.  
2. Special care should also be taken to remove imaging powders very thoroughly with pumice and water.  
3. APPLICATION OF G-CEM ONE ADHESIVE ENHANCING PRIMER (Optional)  
Dispense a few drops into a clean dispensing dish. Close the bottle cap immediately.  
b) IMMEDIATELY apply to the prepared enamel and dentin surfaces using the disposable applicator.  
c) Leave undisturbed for 10 seconds after the end of application.  
d) Dry thoroughly for 5 seconds with oil free air until MAXIMUM air pressure. Use vacuum suction to prevent splatter of the adhesive.  
**RESTORATION PREPARATION**  
Make sure that the restoration is pretreated and handled according to the manufacturer's instructions.  
a) For Glass Ceramics, hybrid ceramics and composite, apply a silane coupling agent to the bonding surface, such as G-Multi PRIMER and dry with an air syringe.  
a) Remove the cartridge cover (Fig. 5)  
b) Depress the lever to dispense the required amounts of paste onto the mixing pad.  
NOTE:  
1. Before dispensing the first time from a new cartridge or after a long interval in between use, bleed the pastes from the cartridge. (Fig. 6)  
2. The adjusting guide can be moved forwards or backwards to change the quantity of material to be dispensed. (Fig. 7)  
3. Position of adjusting guide for one mix  
4. Position of adjusting guide for two mix

c) While depressing the lever (Fig. 8-1), level and cut off the extruded paste from the cartridge tips as shown. Clean the cartridge tip and replace the cartridge cover (Fig. 9)

d) When the lever is released after dispensing, residual pastes at the cartridge tips are retracted, resulting in a clean tip. Clean the cartridge tip and replace the cartridge cover (Fig. 9)

**NOTE:**  
1. The amount of residual material in the cartridge can be confirmed by the position of the sliding block on the underside of the dispenser.  
**5. MIXING**  
After curing, spread the two pastes thinly on the mixing pad using a plastic spatula. Mix thoroughly, with lapping strokes, for 10 seconds. Be careful not to incorporate air bubbles.  
a) Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes 30 seconds until mix at 23°C (73°F). Higher temperatures will shorten the working time.  
**NOTE:**  
1. Ambient light may shorten the working time. It is recommended to diminish the intensity of turn off your operation light during the cementation procedure.  
2. When the cement is in contact with G-CEM ONE ADHESIVE ENHANCING PRIMER, the setting reaction is accelerated.  
3. In case of veneer cementation, it is recommended to cement a maximum of one veneer at a time due to relatively short working time.  
**6. MAINTAIN MODERATE PRESSURE** until final polishing is done.  
**7. EXCESS CEMENT REMOVAL**  
a) Maintain moderate pressure and make sure the restoration remains in place.  
b) Remove excess cement using one of the techniques below:  
1. Tack cure with light curing unit.  
2. Apply vacuum suction to remove the excess cement for 1 second until it reaches a solid rubbery consistency and remove the excess cement using a probe or another appropriate instrument. Reaching a good consistency for excess cement removal depends on the light curing unit being used. It is important to find the best combination of time, intensity, and distance of the curing light in use.  
3. Keep moderate pressure and remove excess cement using a brush. Make sure that excess cement is thoroughly removed from areas of difficult access (e.g., interproximal areas).  
**NOTE:**  
1. If the excess cement is removed before complete hardening, a suitable instrument must be used to hold the restoration in position.  
2. To avoid the oxygen inhibition layer at the margins of the restoration, protect the margins with a glycerin gel after removing the excess cement.  
**8. FINAL SET**  
While maintaining moderate pressure, light cure all surfaces / margins. Refer to the following chart for irradiation time.

Irradiation time	
High power LED (more than 1200 mW/cm <sup>2</sup> )	10 seconds
Halogen / LED (700-1200 mW/cm <sup>2</sup> )	20 seconds

In case of the following restorations, let the material set for 4 minutes after removing excess cement:  
- Metal and metal coping restorations  
- Metal-free dark and opaque restorations  
- Translucent metal-free restorations greater than 2 mm in thickness

**9. FINAL POLISHING**  
If needed, finish and polish margins with appropriate instruments.

**B. CEMENTATION TECHNIQUE for metal, ceramic, fiber posts, and cast post and cores**

1. TOOTH PREPARATION  
Prepare the post space according to manufacturer's instructions. The post space should be cleaned, rinsed and dried thoroughly using paper points. The post space should be dry. If extra adhesion is needed, apply G-CEM ONE ADHESIVE ENHANCING PRIMER to adhesive surface for 10 seconds. Dry thoroughly with oil-free air until MAXIMUM pressure until the primer solution does not come out of the root canal entrance anymore. Use vacuum suction to prevent splatter of the adhesive. Remove excess primer solution with paper points.  
**NOTE:**  
1. Sodium hypochlorite solution is recommended to use for the chemical cleaning of the post space.  
2. Do not use H<sub>2</sub>O<sub>2</sub> and/or EDTA to chemically clean the post space as these products may lower the bond strength.  
**3. POST PREPARATION**  
Make sure that the post is pretreated and handled according to the manufacturer's instructions.  
**4. DISPENSING & MIXING**  
Refer to A. 4 and A. 5.  
**4. CEMENTATION**  
a) Apply the cement on adhesive surface of the post.  
b) Place the post immediately into the post space. Be careful not to incorporate air bubbles enter the root canal.  
**NOTE:**  
1. The oral temperature accelerates the setting reaction of the cement.  
2. Ambient light may shorten the working time. It is recommended to diminish the intensity of turn off your operation light during the cementation procedure.  
3. DO NOT use lentulo spirals to place the cement in the post space as this may shorten the working time.  
4. When the cement is in contact with G-CEM ONE ADHESIVE ENHANCING PRIMER, the setting reaction is accelerated. In order to make excess cement removal easier, make sure proper air drying and not to lean on excess portion of G-CEM ONE ADHESIVE ENHANCING PRIMER around adhesive area.  
c) Maintain moderate pressure.  
**5. EXCESS CEMENT REMOVAL**  
Continue to maintain moderate pressure making sure the post remains in place, and remove excess cement with appropriate instruments. In case of a cast post and core, tack cure all surface for 1 second each with a visible light curing unit.  
**6. FINAL SET**  
a) While maintaining moderate pressure, light cure all surfaces / margins. Refer to the following chart for irradiation time.  

Irradiation time	
High power LED (more than 1200 mW/cm <sup>2</sup> )	10 seconds
Halogen / LED (700-1200 mW/cm <sup>2</sup> )	20 seconds

  
b) Let the material set for another 4 minutes.

**CLEANING / STERILIZATION OF PASTE PAK DISPENSER**  
For cleaning and sterilizing of PASTE PAK DISPENSER, refer to the respective instructions for use.

**STORAGE**  
Recommended for optimal performance, store at temperature of 4-25°C (39.2-77.0°F) away from heat, moisture and direct sunlight.

**SHADE**  
A2 (Vita® shade), A03 (opaque), White Opaque, Translucent  
Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

**PACKAGES**  
G-CEM ONE Paste Pak Starter kit (A2, Translucent);  
G-CEM ONE Paste Pak Cartridge (7.6 mL) (1); G-CEM ONE ADHESIVE ENHANCING PRIMER (2 mL) (1); Mixing Pad (No.22) (1)  
G-CEM ONE Paste Pak Single Refill (A2, A03, White Opaque, Translucent);  
G-CEM ONE Paste Pak Cartridge (7.6 mL) (1); Mixing Pad (No.22) (1)

**CAUTION**  
1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton swab. Flush with water. To avoid contact, a rubber dam and/or COCCA BUTTER can be used to isolate the operation field from the oral or tissue surface.  
2. In case of contact with eyes, flush immediately with water and seek medical attention.  
3. Do not use G-CEM ONE Paste Pak in combination with eugenol containing materials as eugenol may hinder G-CEM ONE Paste Pak curing or bonding.  
4. Use a protective light shield or personal protective eye wear during light curing.  
5. This product is not indicated for filling or core build-up.  
6. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.  
7. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.  
8. Avoid inhalation or ingestion of material.  
9. G-CEM ONE ADHESIVE ENHANCING PRIMER is volatile. Use in a well-ventilated place.  
10. Do not mix with other products.  
11. Light cure completely for an effective bonding strength. Lower light intensity may cause insufficient adhesion.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://sea.gcsasidental.com>

They can also be obtained from your supplier.  
**CLEANING AND DISINFECTING**  
MULTI-USE DELIVERY SYSTEMS: To avoid cross-contamination between patients this device requires immediate disinfection immediately after use inspect device and label for deterioration. Discard device if damaged. DO NOT IMMERGE: Thoroughly clean device to prevent drying and accumulation of contamination. Disinfect with a mid-level registered healthcare-grade infection control product according to regional / national guidelines.

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