## GC Fuji IX GP FAST

RADIOPAQUE POSTERIOR GLASS IONOMER RESTORATIVE CEMENT IN CAPSULES

GC Fuji IX GP FAST is a fast-setting version of GC Fuji IX GP CAPSULE.

For use only by a dental professional in the recommended indications.

### RECOMMENDED INDICATIONS

- Class I and II restorations in deciduous teeth.
- Non-load bearing Class I and Class II restorations in permanent teeth.
- Intermediate restorative and base material for heavy stress situation in Class I and Class II cavities using sandwich laminate technique.
- Class V and root surface restorations.
- Core build-up.

### CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE	Fast	Original
Powder / Liquid Ratio (g/g)	0.40/0.11	0.40/0.11
Mixing Time (sec.)*	10"	10"
Time of initial extrusion after mixing	10 seconds maximum	
Working Time (min., sec.)*	1'15"	2'00"
Net Setting Time (min., sec.)*	2'00"	2'20"
Final Finishing Commencing Time (37°C, after starting mix) (min., sec.)	3'00"	6'00"

Test conditions : Temperature (23 +/-1°C),

Relative humidity (50 +/-5%)

\*ISO 9917-1 : 2003 (E) (Glass polyalkenoate cement)

### 1. CAVITY PREPARATION

- Prepare tooth using standard techniques. Extensive mechanical retention is unnecessary. For pulp capping use calcium hydroxide.
- Apply GC CAVITY CONDITIONER (10 seconds) or GC DENTIN CONDITIONER (20 seconds) to the bonding surfaces using a cotton pellet or sponge (Fig. 1).

- Rinse thoroughly with water. Blot away excess water with a cotton pellet or dry by gently blowing with an air syringe. DO NOT DESICCATE. Prepared surfaces should appear moist (glistening).

### 2. CAPSULE ACTIVATION AND MIXING

- Before activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 2).
- To activate the capsule, push the plunger until it is flush with the main body (Fig. 3).
- Immediately place the capsule into a metal GC Capsule Applier and click the lever once (Fig. 4). The capsule is now activated.

Note :  
The capsule should be activated just before mixing and used immediately.

- Immediately remove the capsule and set it into a mixer (or an amalgamator) and mix for 10 seconds (+/- 4.000RPM) (Fig. 5).

### 3. RESTORATIVE TECHNIQUE

- Immediately remove the mixed capsule from the mixer and load it into the GC Capsule Applier.
- Make two clicks to prime the capsule then syringe (Fig. 6). The working time is 1 minute 15 seconds from start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.
- Within 10 seconds maximum after mixing, start to extrude the mixture directly into the preparation (Fig. 7).
- Form the preliminary contour, and cover with a matrix if required.
- During the first 3 minutes from start of mix extra care should be taken to avoid moisture contamination or drying-out. In case this cannot be guaranteed, immediately apply or GC Fuji Varnish (blow dry) or GC Fuji Coat LC (light cure).

Note :

- To adjust the direction of the nozzle, hold the applicator with the capsule towards you and turn the capsule body.
- To remove the used capsule, push the applicator release button. Twist the capsule and pull upwards.

### 4. FINISHING

- Final finishing under water spray using standard techniques can begin approximately 3 minutes after start of mixing (Fig. 9).
- Apply a final coat of GC G-COAT PLUS and light cure for 20 seconds (Fig. 9,10). Alternatively, apply GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure) Refer to the respective instructions for use.
- Instruct the patient not to apply pressure for 1 hour.

### STORAGE

Store in a cool and dark place (4-25°C) (39.2-77.0°F). (Shelf life : 2 years from date of manufacture).

### SHADES

A1, A2, A3, A3.5, B2, B3, C4

Shade numbers according to Vita® shade guide.

Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

### PACKAGES

- Single shade package : 50 capsules.
- Assortment package : 50 capsules (10 each of A2, A3, A3.5, B2, B3). Average contents per capsule : 0.40g powder and 0.11g (0.09mL) liquid. Minimum net volume of mixed cement per capsule : 0.14mL. Option : GC CAPSULE APPLIER (1 piece).

### CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.

Last revised : 07/2012

**CE 0086**

**GC**

MANUFACTURED by  
GC CORPORATION  
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

DISTRIBUTED by  
GC CORPORATION  
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU: GC EUROPE N.V.  
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,  
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

GC AMERICA INC.  
3727 West 127th Street, Alsip, IL 60803 U.S.A.  
TEL: +1 708 597 4900

GC ASIA DENTAL PTE. LTD.  
19 Loyang Way, #06-27 Singapore 508724  
TEL: +65 6546 7588

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GC Fuji IX GP FAST Prior to use, carefully read the instructions for use.

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### RECOMMENDED INDICATIONS

- Class I and II restorations in deciduous teeth.
- Non-load bearing Class I and Class II restorations in permanent teeth.
- Intermediate restorative and base material for heavy stress situation in Class II restorations using sandwich laminate technique.
- Class V and root surface restorations.
- Core build-up.

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### DIRECTIONS FOR USE

	Fast	Original
Powder / Liquid Ratio (g/g)	0.40/0.11	0.40/0.11
Mixing Time (sec.)*	10"	10"
Time of initial extrusion after mixing	10 seconds maximum	
Working Time (min., sec.)*	1'15"	2'00"
Net Setting Time (min., sec.)*	2'00"	2'20"
Final Finishing Commencing Time (37°C, after starting mix) (min., sec.)	3'00"	6'00"

Test conditions : Temperature (23 ±1°C), Relative humidity (50 ±5%)

\*ISO 9917-1 : 2003 (E) (Glass polyalkenoate cement)

### 1. CAVITY PREPARATION

- Prepare tooth using standard techniques. Extensive mechanical retention is unnecessary. For pulp capping use calcium hydroxide.
- Apply GC CAVITY CONDITIONER (10 seconds) or GC DENTIN CONDITIONER (20 seconds) to the cavity surfaces using a cotton pellet or sponge (Fig. 1).
- Rinse thoroughly with water. Blot away excess water with a cotton pellet or dry by gently blowing with an air syringe. DO NOT DESICCATE. Prepared surfaces should appear moist (glistening).

### 2. CAPSULE ACTIVATION AND MIXING

- Before activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 2).
  - Push the capsule into the plunger until it is flush with the main body (Fig. 3).
  - Immediately place the capsule into a metal GC Capsule Applier and click the lever once (Fig. 4). The capsule is now activated.
- Note :  
The capsule should be activated just before mixing and used immediately.

### 3. RESTORATIVE TECHNIQUE

- Immediately remove the capsule from the mixer (or an amalgamator) and mix for 10 seconds (+/- 4.000RPM) (Fig. 5).
- Make two clicks to prime the capsule then syringe (Fig. 6). The working time is 1 minute 15 seconds from start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.
- Within 10 seconds maximum after mixing, start to extrude the mixture directly into the preparation (Fig. 7).
- Form the preliminary contour, and cover with a matrix if required.
- During the first 3 minutes from start of mix extra care should be taken to avoid moisture contamination or drying-out. In case this cannot be guaranteed, immediately apply GC Fuji Varnish (blow dry) or GC Fuji Coat LC (light cure).

### 4. FINISHING

- Final finishing under water spray using standard techniques can begin approximately 3 minutes after start of mixing (Fig. 9).
- Apply a final coat of G-COAT PLUS and light cure for 20 seconds (Fig. 9, 10).  
Alternatively, apply GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure) Refer to the respective instructions for use.
- Inform the patient not to apply pressure for 1 hour.

### STORAGE

Store in a cool and dark place (4-25°C) (39,2-77,0°F). Shelf life : 2 years from date of manufacture.

### SHADES

A1, A2, A3, A3.5, B2, B3, C4  
Shade numbers according to Vita® shade guide.  
Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

### PACKAGES

- Single shade package : 50 capsules.
- Assortment package : 50 capsules (10 each of A2, A3, A3.5, B2, B3).  
Average content per capsule : 0.40g powder and 0.11g (0.09mL) liquid.  
Minimum net volume of mixed cement per capsule : 0.14mL.  
Option : GC CAPSULE APPLIER (1 piece).

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