## Saliva-Check BUFFER

In Vitro Test for checking the Quality, pH and Buffering Capacity of Saliva

For use only by a dental professional in the recommended indications.

### RECOMMENDED INDICATIONS

To check the properties of resting and stimulated saliva.

### CONTRAINDICATIONS

In vivo diagnostic use.

### DIRECTIONS FOR USE

Prior to any visit where a saliva diagnostic procedure is planned, instruct the patient not to smoke, consume food or drink, brush the teeth or use a mouth wash for at least one hour prior to the scheduled appointment time.

#### I. TESTING OF RESTING SALIVA

**TEST 1 – Visual inspection of level of hydration**  
Visually assess the lower lip labial gland secretion. Evert the lower lip, gently blot the labial mucosa with a small piece of gauze and observe the mucosa under good light. Droplets of saliva will form at the orifices of the minor glands.

Assess the time for visible production of saliva as follows<sup>1</sup>:

Greater than 60 seconds : resting flow **Low** ●  
Less than 60 seconds : resting flow **Normal** ●


#### TEST 2 – Saliva consistency

Visually assess the resting salivary consistency in the oral cavity<sup>1</sup>.

Sticky frothy saliva residues: Increased viscosity ●  
Frothy bubbly saliva: Increased viscosity ●  
Watery clear saliva: Normal viscosity ●

#### TEST 3 – pH measurement

Instruct the patient to expectorate any pooled saliva into the collection cup. Take a pH test strip, place this into the sample of resting saliva for 10 seconds, and then check the colour of the strip. This should be compared with the testing chart available in the package.

#### II. TESTING OF STIMULATED SALIVA

##### TEST 4 – Saliva quantity

Instruct the patient to chew the piece of wax to stimulate salivary flow. After 30 seconds, let the patient expectorate into the spittoon. Continue chewing for a further 5 minutes, collecting all the saliva into the collection cup at regular intervals.

The quantity of saliva can be measured by checking the mL markings on the side of the cup.

Quantity of saliva at 5 minutes<sup>2</sup>

< 3.5 mL **Very low** ●  
Between 5.0 – 3.5 mL **Low** ●  
> 5.0 mL **Normal** ●


Note: Normal stimulated saliva flow rate may vary between 1mL/min – 1.6mL/min.

##### TEST 5 – Buffering capacity

- Remove a Buffer test strip from the foil package and place onto an absorbent tissue with the test side up.
- Using a pipette, draw sufficient saliva from the collection cup and dispense one drop onto each of the 3 test pads. Immediately turn the strip 90° to soak up excess saliva on the absorbent tissue. This will prevent the excess saliva from swelling on the test pad and possibly affecting the accuracy of the test result.
- The test pads will begin to change colour immediately and after 2 minutes the final result can be calculated by adding the points according to the final colour of each pad. See conversion table and examples underneath.

Conversion table	
Test pad colour at 2 minutes	
Green	4 points
Green/Blue	3 points*
Blue	2 points
Red/Blue	1 point *
Red	0 points

#### Examples

\*Where a colour combination provides an unclear result, use intermediate scores.

#### Interpreting the result.

Combined total **Buffering ability of Saliva**

0-5 **Very low** ●  
6-9 **Low** ●  
10-12 **Normal / High** ●

### STORAGE

Store components at room temperatures (15° – 30°C). (Shelf life: 3 years from date of manufacture)

### PACKAGE

Contents	Quantity
In vitro pH test strips	20
Saliva collection cup	20
Wax for saliva stimulation	20
Saliva dispensing pipette	20
Buffer test strip	20

### CAUTION

- In warm ambient temperatures, store the wax in a refrigerator.
- After opening the wax and Buffer test strip, use immediately.
- All components used for each test are for single use only. Do not reuse.
- Dispose of according to local regulations.

<sup>1</sup> Information by kind permission of Prof. L Walsh©2002

<sup>2</sup> Lund University Faculty of Odontology Department of Cariology, Sweden

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